



## Safety Storage Cabinets

- Flammable Storage Cabinet (Yellow)
- Combustible Storage Cabinet (Red)
- Corrosive Storage Cabinet (Blue)
- Under counter / Fume Hood Storage Cabinet
- Designed in accordance to US National Fire Protection Agency Code 30

## Special Features

- Constructed from industrial grade electrogalvanised steel 18 gauge for strength and durability
- Epoxy powder coated to meet harsh industrial / laboratory environments
- Double layered construction with 1.5" airspace between layers
- Fully welded 2" sill to contain chemical spills within the cabinet
- 2 ventilation ports protected with flash arrestors
- 3 point door latching mechanism
- Heavy duty adjustable foot pads for leveling
- Large Warning Labels for easy identification
- Adjustable shelves with 2" centers



# Specifications

MODELS (gallon)	15	30	45	60	90
<b>Dimensions L x W x H (inches)</b>	24 x 21.625 x 35.75	43 x 18 x 44	43 x 18 x 65	34 x 34 x 65	43 x 34 x 65
<b>Weight</b>	75 Kg	105 Kg	145 Kg	180 Kg	220 Kg
<b>Design</b>	Manual closing door				
<b>Door</b>	1 door	2 door	2 door	2 door	2 door
<b>Type / Colour</b>	Flammable – Yellow Combustible – Red Corrosive – Blue Hazardous – Grey				
<b>Shelves</b>	1	2	2	2	2
<b>PP Lining</b>	Only for Corrosive (Optional for the rest)				

Ordering: please quote Model No, followed by type and PP if lining required.



**AdvanceLab (S) Pte Ltd**

22 Senang Crescent, Singapore 416595

Tel: +65 6448 8255 Fax: +65 6448 9833

[www.advancelab.com.sg](http://www.advancelab.com.sg)

Email: [info@advancelab.com.sg](mailto:info@advancelab.com.sg)

